

ABSORBANCE MEASUREMENT ON μ Drop PLATE
Introduction

The Thermo Scientific™ μ Drop™ Plate is a quick and easy tool for DNA and RNA assays using photometric measurement and sample volumes down to 2 μ L.

The μ Drop Plate provides a straightforward way of analyzing up to 16 microliter-scale samples simultaneously. The fixed light path of the μ Drop Plate allows direct calculation of the nucleic acid concentrations of the samples. Samples are easy to pipette onto the plate with a single- or an eight-channel pipette. It takes less than 10 seconds to measure the full DNA spectra from all 16 samples with **Hidex Sense** platereader. For flexibility, the μ Drop Plate also enables cuvette measurement. The μ Drop Plate is available from Hidex.

The **Hidex Sense** will help your lab become more effective. The touch screen user interface makes the operation safe and comfortable. Straightforward application focused operation minimizes time spent on instrument training, and is essential for superior results. ELISA, protein quantification and enzyme activity assays are measured with ultrafast full spectrum readout, using

Materials *DNA* Calf Thymus DNA, Lyophilized Powder, Type XV (Sigma D 4522)
 Microplate μ Drop (Hidex Code 425-5801)


Hidex Sense

 μ Drop Plate
Test procedure


5 μ L of DNA dilutions were pipetted in duplicates into the measurement positions on the plate. Absorbances were measured at 260 nm and at 320 nm. The 320 nm OD values were used as blanks.

Reader setup

Flashes 10
 Wavelengths Center 260 nm (window 5)
 Center 320 nm (window 5)

Results

Conc (mg/L)	260 nm			320 nm (Blank)			Avg. - Blank
	Repl 1	Repl 2	Avg.	Repl 1	Repl 2	Avg.	
50	0.129	0.125	0.127	0.043	0.044	0.044	0.084
125	0.245	0.215	0.230	0.038	0.038	0.038	0.192
250	0.361	0.347	0.354	0.042	0.049	0.046	0.309
500	0.658	0.651	0.655	0.043	0.040	0.042	0.613
1000	1.298	1.285	1.292	0.046	0.042	0.044	1.248



Absorbance vs. DNA concentration.

Conclusions

μ Drop Plate is a practical tool in photometric DNA and RNA quantitation and purity analysis on Hidex Sense platereader. The spectrograph used in absorbance measurement enables simultaneous detection several wavelengths. The DNA purity test is a good example about multi-wavelength assay in which OD values at 260 nm, 280 nm and 320 nm are required.