

# Evaluation of Active Counting Mode option in LabLogic

## $\beta$ -RAM Model 5 and Laura software


J Attwood, B O'Brien (AstraZeneca R&D Charnwood); H W Loaring, A Williams (LabLogic Systems Ltd.)

### Introduction



LabLogic have recently introduced the  $\beta$ -RAM<sup>®</sup> Model 5 radio-HPLC detection system, which includes the option of Active Counting Mode (ACM<sup>™</sup>), a new online counting system with “unrivalled limits of detection” and “improved resolution and peak definition” compared to older systems. ACM enhances counting conditions in real time. Following the detection of a peak, ACM optimises scintillant flow rate, which impacts upon analyte residence time and improves signal-to-noise. We have evaluated this new technology by assessing the comparability of online detection between standard radio-HPLC and ACM radio-HPLC, as well as standard radio-UPLC and ACM radio-UPLC.



### ACM Parameters



----- Background subtraction (calculated pre-injection)




### Method


A pooled human urine sample with nominal <sup>14</sup>C concentrations of 18,000 DPM/75 $\mu$ L (HPLC) and 2,400 DPM/10 $\mu$ L (UPLC) was utilised as a source of radioactivity for all analyses. Sub-aliquots of this sample were concentrated by freeze-drying or diluted where appropriate, prior to Liquid Scintillation Counting (LSC) to confirm <sup>14</sup>C content. HPLC: Waters Xterra MS C18 column; UPLC: Acquity UPLC<sup>®</sup> C18 column; HPLC & UPLC: gradient elution using ammonium bicarbonate and acetonitrile:methanol with 2:1 scintillant:eluate flow rate. The following parameters and settings were assessed: radio-HPLC & radio-UPLC: Background Threshold (0,1,2,3,5 & 10), Scintillant Flow Multiplier (x2 and x3) and Flush Time (4, 8 & 16 seconds); radio-UPLC only: flow-cells (200, 500 & 1000  $\mu$ L); radio-HPLC: 500  $\mu$ L flow-cell only.

### Results

#### HPLC-ACM: Background Threshold




#### HPLC-ACM: Flush Time




#### UPLC-ACM: Flow Cells

As expected, increasing flow-cell size resulted in broadened peaks and decreased peak resolution due to a lower residence time of the analyte

#### HPLC-ACM: Scintillant Flow Multiplier



#### HPLC-ACM vs. UPLC-ACM



### Discussion

- Reproducibility was generally good, although detection of minor peaks varied between some replicates.
- Flush time increased peak height but may merge neighbouring peaks into one peak.
- UPLC-ACM requires less sample and displays better limits of detection than HPLC-ACM.
- No observable difference between Background Threshold 2 and above.
- Scintillant flow rate setting should be optimised to the methodology.
- Further investigations into the comparability of ACM and off-line radioactivity counting, as well as ACM analysis of <sup>3</sup>H are required.

### Conclusions

- ACM should be used with optimised chromatography.
- Peak reject criteria needed for smaller, insignificant peaks.
- Background threshold could be scaled down from 0-10 to 0-2.