

HIDEX



Hidex 300 SL

Liquid Scintillation Counter

www.lablogic.com

 **LabLogic**
EXPERIENCE & EXPERTISE

The most advanced transportable and user friendly LSC on the market

The Hidex 300 SL is a revolutionary instrument which incorporates the most advanced technology available, triple-to-double coincidence ratio (TDCR) counting.

Compact and Transportable



The 300 SL has a modern and compact design measuring half the size and weight of some of its long standing rivals. It is therefore much easier to fit into smaller, more space conscious laboratories, such as transporting in research vessels or mobile labs.

User Friendly Software


The 300 SL uses an external Windows 10 PC, enabling several users to transfer their method and data files to other offline or networked PC's.

The MikroWin 300 SL software features an easy-to-use graphical interface and advanced data reduction capabilities. An unlimited number of protocols, easy export to Excel™, custom calculations and quench curves are all included.

The 300 SL can also operate with Laura software for 21 CFR Part 11 compliance for secure handling of method and data files.



Sample ID	Activity	Parameter 1	Parameter 2	Parameter 3	Parameter 4	Parameter 5	Parameter 6	Parameter 7	Parameter 8
1	100	100	100	100	100	100	100	100	100
2	200	200	200	200	200	200	200	200	200
3	300	300	300	300	300	300	300	300	300
4	400	400	400	400	400	400	400	400	400
5	500	500	500	500	500	500	500	500	500
6	600	600	600	600	600	600	600	600	600
7	700	700	700	700	700	700	700	700	700
8	800	800	800	800	800	800	800	800	800



TDCR Technology


The 300 SL utilises a unique design with three photomultiplier detectors aligned 120 degrees from each other, which offers superior efficiency and counting results for both advanced research and environmental work.

In addition, the 300 SL can be used in triple coincidence mode only, which removes interference from chemiluminescence and the need to dark adapt. The counter can also be operated in dual coincidence mode using conventional standards and quench correction methods.

Quench curves can be generated and stored within the MikroWin software, which can be automatically applied or manually fitted to recalculate previous data.

No internal or external standard source

The TDCR method allows easy determination of sample counting efficiency without the use of any internal or external standard source of radioactivity, offering many benefits for the modern lab including no hidden disposal costs, no health hazard when servicing, and no transportation restrictions.



Advanced spectrum analysis highlighting triple and double count spectrums for each sample.

Optimal Lead Shield Design

A minimum of 70 mm lead in all directions provides excellent shielding and minimises instrument weight.

Copper Shielding

Eliminates X-rays from lead shield.

Measurement Chamber

Highly reflective proprietary paint maximises photon detection.

Lead Shutter

Reduces cosmic radiation.

Robotic Loading Arm


Removes the need for a complex elevator mechanism, preventing jams, and is easily serviceable.

No Source Required

Eliminating hazards and cost without compromising results.

PMTs

Three PMTs detectors in coincidence mode provide optimal counting geometry and facilitate TDCR calculations.



Super Low-Level Model

Similar to the standard instrument the Hidex 300 SL super low-level scintillation counter is equipped with additional lead shielding, low-level PMT detectors, and an active guard detector for further background reduction.


Ideal for ³H in water measurements, other low-level environmental monitoring, radiocarbon dating and biofuel verification.

New Active Guard

The Active Guard is a separate scintillator detector which measures and subtracts real-time background radiation. The guard can be turned off for high energy samples and does not interfere with alpha/beta separation.

300 SL Options

- Alpha/Beta Separation Module
- Temperature Control Module
- Low-Level PMT Detectors
- Digital Lead Shield
- External Eu-152 Standard



2D/3D spectral separation of alpha's from beta's based on pulse length discrimination.

Specifications	
Size	63 (L) x 52 (W) x 68.5 (H) cm
Weight	Standard Model: 130 kg Super Low Level Model: 180 kg
Electrical Connections	100 - 240 V 50 - 60 Hz
Energy Range β's	0 - 2,000 keV
Energy Range α's	0 - 10,000 keV
Efficiency	H-3 unquenched 68 - 72% H-3 (8 mL water, 12 mL cocktail) > 30% C-14 unquenched 96 - 98% α's (Po-210, U-234/U238, Am-241, Rn-222, Ra-226) >95%

Please refer to the Technical Specification Sheet for further information

600 SL Automatic TDCR Liquid Scintillation Counter

Designed to meet the needs of laboratories processing large quantities of samples, the new Hidex 600 SL is a high throughput automatic TDCR liquid scintillation counter.

Proven technology

The Hidex 600 SL uses the robust and unique triple-to-double coincidence ratio (TDCR) counting technology from the successful 300 SL series. Coupled with added sample capacity for up to 500 small vials (or 210 large vials), the 600 SL can process samples at a rate which will satisfy even the most demanding laboratories.

The instrument's software allows the user to work with an unlimited number of method files. Methods are barcode identified for automatic processing. Data files can be automatically exported in Excel®, csv or text for integration into LIM systems.

Integration with Laura software

The instrument can also be linked to Laura, the industry-standard radiochromatography software, which provides all the tools to enable the counter to be used in a 21 CFR part 11 compliant manner.

A major benefit is that time-stamped audit trails are generated automatically and cannot be modified, making it easier to track end-user activity.

Also in line with the data integrity guidelines, each end-user will have their own unique, secure ID login and password. Different levels of access can also be granted to specific end users.

Database data storage option

Allows database only storage of precious data. Laura supports both Oracle and MSSQL, with the data only accessible through Laura.

Multi-level security

Definable levels of access are controlled for each user.

Project level access

User access can be defined on a project basis.

Configurable audit trail


A comprehensive audit trail provides full accountability for data and user actions.

Electronic signatures

Configurable e-signatures are available, including dual signatures if required.

Lock projects, methods, evaluations against editing

Areas can be secured against editing using the user's login information.



Optional Features

The Hidex 600 SL is available with all the options of the standard model, such as powerful alpha/beta separation, low level PMT detectors, cooling unit, and internal Eu-152 reference source.

Specifications	
Size	125 (W) x 69 (H) x 64 (D) cm
Weight	200 kg
Sample Capacity (20 mL/7 mL)	210/500
Background (CPM)	12 (H-3 in water)
Counting Efficiency (%)	H-3 >70% C-14 >96%

Please refer to the Technical Specification Sheet for further information

Service and Support

Users of our systems can benefit from our comprehensive, fully inclusive service and support.

We can give reassurance that if things go wrong or you need expert advice, help is only an e-mail or phone call away.

Validation Services

Our Validation Service enables you to implement and get maximum value from your investments as soon as possible.

We work as a partner with your Quality Manager, System Manager and users to provide a tailored Validation Plan, suited to your needs. Our Validation Specialists have years of experience in GLP system validation, detailed knowledge of our systems, together with other industry standard systems to help you meet company and regulatory requirements.

Training

LabLogic can provide a variety of training courses and workshops to help you get the most out of your instrument and software.

All training is performed by our expert Product and Support Specialists who have many years experience in the development and use of the instruments and software.

Certificates can be provided to complement your internal GLP training records.

Europe & Worldwide

LabLogic Systems Limited

Paradigm House, 3 Melbourne Avenue
Broomhill, Sheffield, S10 2QJ, UK

E-mail: solutions@lablogic.com

Tel: +44 (0)114 266 7267

Fax: +44 (0)114 266 3944

www.lablogic.com

USA & Canada

LabLogic Systems, Inc.

1911 N US HWY 301, Suite 140
Tampa, FL 33619, USA

E-mail: solutions@lablogic.com

Tel: +1-813-626-6848

Fax: +1-813-620-3708

www.lablogic.com



LabLogic

EXPERIENCE & EXPERTISE



Certificate No: 1535
ISO 9001



UKAS
MANAGEMENT
0526

Silver
Until 2022



Certificate No: 10926
ISO 9001