# Flow-RAM\*

PET/SPECT radio-HPLC Detector





The innovative, flexible and compliant flow-through detector

LabLogic's PET/SPECT radio-HPLC Detector, the innovative Flow-RAM™, is designed to meet the ever increasing demands of the modern nuclear medicine laboratory.

# **Flexibility**

The Flow-RAM™ is available in a single or dual detector configuration. In the latter configuration, each of the detectors can be set up independent of each other and used for different applications, if required.

With a wide range of detectors suitable for both low and high levels of radioactivity, the Flow-RAM™ is an ideal quality control and research instrument.

# Compliance

As with all LabLogic products, the Flow-RAM™ has been designed to exceed regulatory requirements. The instrument's front panel only has a power button and display so all detector parameters are controlled via the software method. Used with our industry standard Laura for PET™ software, the Flow-RAM™ ensures regulatory compliance via controlled access and audit-trail.

# Controlled by Laura for PET™

Flow-RAM<sup>™</sup> control, digital data collection, analysis and reporting is provided by Laura for PET™, the industry standard radiochromatography data system developed by LabLogic.



# NaI Well-type Detector

### A Range of Detectors

Suitable for PET, SPECT and High Energy Beta radionuclides across a range of radioactivity levels.

# You felaton Heavenent Jos Minton Heb

# **USB** Interface

LabLogic

**Single Power Button** 

Flow-RAM\*

Operating parameters such as high voltage, upper and lower level discriminators as well as power supply and data transfer are all achieved through a USB connection.

# **Built-in Analog to Digital Converter**

This feature coverts analog signals from other detectors such as UV, ECD, etc., into digital for use within Laura™, bringing all signals into one place.

### **Built-in Analog Output**

Users can take an analog signal into any third party software.

### **Bespoke Lead Shielding**

In order to minimize background interference LabLogic provides a range of lead shielding to meet different requirements.

### **Intelligent Front End Display** Connectivity

The Flow-RAM™ is compatible with all commercially available HPLC systems.

EXPERIENCE & EXPERTISE

Users are able to view instrument parameters and performance at a glance.

# Basic Specifications

7.5 (H) x 22 (L) x 23 (W) cm Size Weight 4.5 lbs Connectivity USB

Please refer to the Technical Specification Sheet for further information

Flow-RAM® EXPERIENCE & EXPERTISE

# Dual Flow-RAM™

The Dual Flow-RAM™ can have up to two detectors connected to it which can be used simultaneously and independently, thus providing maximum flexibility and a two-in-one solution.

# A Range of Detectors

Suitable for PET, SPECT and High Energy Beta radionuclides across a range of radioactivity levels.



# **Detector Options**

radio-HPLC Det	ector Options	
Detector	Radioactivity Type	Commonly Used Isotopes
1" NaI PMT	PET SPECT	F-18, C-11, Gα-68, Rb-82 Tc-99m, In-111
2" NaI PMT	High Energy Gamma	F-18, C-11, Ga-68, Rb-82, Zr-89
lastic PMT	Beta	Lu-177, Y-90, I-131, Re-188, Re-186
).1" NaI PMT	Low Energy Gamma	I-125
ell-Type NaI PMT	PET or SPECT with low amou i.e. small animal imaging app	nts of activity olications and measuring low-level impurities.
N Diode	Semi-prep HPLC purification	of radio-tracer after synthesis.
sI PIN Diode		and occasionally quality control of sitive than standard PIN Diode.
/ww.lablogic.co	om ———	



# Flow-RAM<sup>™</sup> is controlled by the industry standard radiochromatography software

Flow-RAM™ control, digital data collection, analysis and reporting is provided by Laura for PET™, the industry standard radiochromatography data system developed by LabLogic.

# Single Software Solution

Laura for PET™ is a single software solution for the PET/SPECT QC environment. Rather than having to use several software systems, the QC analyst only needs Laura for PET with radio-HPLC, HPLC, radio-TLC, GC and MCA analysis.

# **Regulatory Compliance**

Laura for PET™ is designed to meet the regulatory compliance requirements for GMP and FDA 21 CFR Part 11 and FDA 21 CFR Part 212 requirements.

### Intuitive

Laura for PET™ is easy to use allowing quick familiarisation with the system and rapid adoption by its users.

#### **Prudent Investment**

From a standalone installation to a full client/server implementation, Laura for PET™ offers a scalable solution to your data storage requirements. Applications and data can be separated across local clients and central servers to match IT requirements.

### **Complete Workflow Solution**

Laura for PET™ will guide the analyst through the main QC process. Methods can be pre-defined and chosen from a drop down list. Automatic settings provide for automatic peak integration and report printing.

### **Consistent and Secure**

The software allows the user to lock the configuration of methods, reports, etc., to ensure consistency of data capture and presentation.



# **Digital Signal**

Signals from radiochromatography instruments are reported digitally in counts, CPS and CPM. In addition, the full dynamic range of the detector is processed avoiding limitations commonly seen with analog signals.

# **Half Life Correction**

Half life correction function enables users to correct for half life while a run is in progress or post-run, using a reference time and date.



### Reporting

Relevant values are reported such as % total for radiochemical purity checks. The powerful Report Designer allows you to configure multiple reports and to calculate and present data in real time.



# Service and Support

Users of our systems can benefit from our comprehensive, fully inclusive service and support.

We can give reassurance that if things go wrong or you need expert advice, help is only an e-mail or phone call away.



# **Validation Services**

Our Validation Service enables you to implement and get maximum value from your investments as soon as possible.

We work as a partner with your Quality Manager, System Manager and users to provide a tailored Validation Plan, suited to your needs. Our Validation Specialists have years of experience in GLP system validation, detailed knowledge of our systems, together with other industry standard systems to help you meet company and regulatory requirements.



# **Training**

LabLogic can provide a variety of training courses and workshops to help you get the most out of your instrument and software.

All training is performed by our expert Product and Support Specialists who have many years experience in the development and use of the instruments and software.

Certificates can be provided to complement your internal GLP training records.



















# USA & Canada LabLogic Systems, Inc.

East Pointe Park, 1040 East Brandon Blvd. Brandon, FL 33511-5509, USA

E-mail: solutions@lablogic.com Tel: +1-813-626-6848

Fax: +1-813-620-3708 Web: www.lablogic.com



# Europe & Worldwide LabLogic Systems Limited

Paradigm House, 3 Melbourne Avenue Broomhill, Sheffield, S10 2QJ, UK

E-mail: solutions@lablogic.com Tel: +44 (0)114 266 7267 Fax: +44 (0)114 266 3944 Web: www.lablogic.com





