

# Hidex Triathler

Liquid Scintillation Counter and Gamma Counter



**HIDEX**

**LabLogic**

EXPERIENCE & EXPERTISE

## Instant results whenever and wherever with this compact and portable counting system

The Hidex Triathler is a single sample counter which provides fast and accurate results for several life science and environmental applications. It can count all radioisotopes including tritium in a variety of sample formats.

Due to its very small size and light weight, Triathler can be taken into the field to measure samples on the spot. Although small, Triathler has many advanced features such as spectrum analysis using a multi-channel analyser (MCA), instant DPM results, single-photon luminescence counting, and optional PC control and data transfer.

### Easy to Use

A keypad allows single key operation for immediate results using preset protocols for any isotope.

### Flexible

The Triathler is suitable for just about any beta radiation, gamma radiation or luminescence application and accepts most types of vials and sample formats.



### Portable

The Triathler's small size and light weight make it ideal for personal use on a benchtop or for on-site field/sea measurements.

### Sophisticated

Triathler has advanced features such as a built-in multi-channel analyser, optional alpha-beta separation electronics, and connectivity to a PC for instrument control and data export.



## Triathler in Monitoring

### Wipe Tests

Triathler provides fast and reliable results for regulatory wipe tests in lab areas.

### Water Measurements

Triathler has alpha-beta separation capabilities which makes it ideal to detect alpha-isotopes like Radon ( $^{222}\text{Rn}$ ) in water. It is also sensitive enough to perform monitoring of  $^3\text{H}$  in discharge water used at power plant sites.

### Soil Measurements

Used in conjunction with recognised extraction techniques, Triathler can measure soil samples for contamination with any radioisotope.



## Triathler in Life Sciences

### Molecular Biology

Triathler is well suited for metabolic studies, genetic studies, cell proliferation assays, receptor-ligand assays, and other applications in biosciences. For example, using the optional plastic scintillator adapter (PSA),  $^{32}\text{P}$  can be measured non-destructively (without cocktail).

### Marine Biology

Triathler is commonly used to measure growth of phytoplankton in seawater and lakes by measuring uptake of  $^{14}\text{C}$ .


### Emergency Response

Triathler is well suited for in-field emergency response, e.g. detecting rapid alpha measurements such as  $^{241}\text{Am}$  as well as  $^{137}\text{Cs}$  and  $^{131}\text{I}$ . The external NaI detector can accommodate up to 1 L Marinelli beakers\* to measure food, soil, and liquids (water, milk, etc.) for gamma contamination.

\* Only available with the Becquerel Finder model of Triathler.

### Research

Triathler is ideal as an easy-to-use diagnostic tool in research for detection of beta and gamma isotopes, such as  $^{125}\text{I}$ ,  $^{18}\text{F}$ ,  $^{137}\text{Cs}$ ,  $^{57}\text{Co}$ , etc.



**Basic Specifications**

Size	33 (L) x 25 (W) x 19 (H) cm
Weight	9 kg
Power	110 - 240 V AC, 12 VDC
Sample Types	LSC Vials, Microtubes, Test Tubes
Detector	Single-Photon Counting PMT
Display	2 x 16 character alpha-numeric LCD
Energy Range	2 keV - 2,000 keV
Counting Time	0.1 seconds - 99999 minutes
Output	RS-232C/USB to PC

Please refer to the Technical Specification Sheet for further information

**Specifications – Liquid Scintillation Counting**

Sample Types	LSC Vials, Microtubes, Test Tubes
Maximum Count Rate	2,000,000 CPM
Beta Efficiency	Up to 45% for $^3\text{H}$
Direct $^{32}\text{P}$ Efficiency	Up to 75% in PSA Tube

**Specifications – Gamma Counting**

Sample Types	Tubes or Vials (up to 13 mm diameter)
Detector	32 x 32 mm NaI (TI) crystal (through-hole)
Background Shield	10 mm lead

**Specifications – External Gamma Counting**

Sample Types	LSC Vials, Microtubes, Test Tubes, Marinelli Beakers*
Detector	50 x 50 mm NaI (TI) crystal (planar type)
Background Shield	30 or 50 mm lead

\*Only available with the Becquerel Finder model of Triathler





## Service and Support

Users of our systems can benefit from our comprehensive, fully inclusive service and support.

We provide complete service and support for all of our customers to give reassurance that if things go wrong or you need expert advice, help is only an e-mail or phone call away.



## Validation Services

Our Validation Service enables you to implement and get maximum value from your investments as soon as possible.

We work as a partner with your Quality Manager, System Manager and users to provide a tailored Validation Plan suited to your needs. Our Validation Specialists incorporate years of experience in GLP system validation, detailed knowledge of our systems, together with other industry standard systems to help you meet your company's requirements.



## Training

LabLogic can provide a variety of training courses and workshops to help you get the most out of your instrument and software.

All training is performed by our expert Product and Support Specialists who have many years experience in the development and use of the instruments and software.

Certificates can be provided to compliment your internal GLP training records.

## Related Products

### Hidex 300 SL

Liquid Scintillation Counter



### Hidex Sense

Multi-technology Microplate Reader




### Hidex AMG

Automatic Gamma Counter



### Wilma

On-line Water Monitor  
for Radionuclide Detection



### Beta-RAM

$\beta$  Radio-HPLC Flow Detector



## Consumables



### Europe & Worldwide

#### LabLogic Systems Limited

Paradigm House, 3 Melbourne Avenue  
Broomhill, Sheffield, S10 2QJ, UK

E-mail: [solutions@lablogic.com](mailto:solutions@lablogic.com)

Tel: +44 (0)114 266 7267

Fax: +44 (0)114 266 3944

[www.lablogic.com](http://www.lablogic.com)



Certificate No: 1535  
ISO 9001



Silver  
Until 2022

### USA & Canada

#### LabLogic Systems, Inc.

1911 N US HWY 301, Suite 140  
Tampa, FL 33619, USA

E-mail: [solutions@lablogic.com](mailto:solutions@lablogic.com)

Tel: +1-813-626-6848

Fax: +1-813-620-3708

[www.lablogic.com](http://www.lablogic.com)



Certificate No: 10926  
ISO 9001