

# Bubble-Point™

Filter Integrity Tester



The LabLogic Bubble Point Filter Integrity Tester is a simple to use and compact system which is compatible with most commonly used filters.

[www.lablogic.com](http://www.lablogic.com)

**LabLogic**

EXPERIENCE & EXPERTISE

## Ensuring the integrity of sterilising filters

The integrity testing of sterilising filters is a fundamental requirement of critical process filtration applications in the pharmaceutical and PET industries. The ultimate measure of a sterilising filter's performance is bacterial retention.

The most widely used non-destructive integrity test is the bubble point test. The bubble point test exploits the fact that liquid is held in the pores of a filter by surface tension and capillary forces. The minimum pressure required to force liquid out of the pores is a measure of the pore diameter. The bubble point test detects minor filter defects and out-of-size pores and correlates with the bacteria passage test.

The bubble point test is a sensitive visual technique performed routinely as part of the suite of product release tests in PET radiopharmaceutical production centers.

## Key Features

- Integrated and compact system.
- Compatible with most filters used for critical process filtration applications.
- Luer lock filter connection for ease and speed of connection.
- Simple air line connection mechanism.
- Calibrated reference gauge.
- Easy-clean anodised surfaces.
- Clearly visible water beaker for identification of bubble point.
- Optional leaded acrylic shield for radiation protection.
- Fine pressure adjustment.
- Can be mounted inside a fume cupboard if necessary.
- Option for complete system IQ/OQ validation.

## Specifications

**Bubble Point Test:** 0 - 6895 mbar (0 - 100 psi)

**Maximum Gas Supply Pressure:** 6895 mbar (100 psi)

**Pressure Gauge Working Temperature Range:**  
-40°C to +60°C (-40°F to +140°F)

**Size:** 12" (W) x 10" (D) x 17" (H) (310 mm x 250 mm x 430 mm)

**Weight:** 7.7 lb (3.5 kg) (without shield), 13 lb (6 kg) (with shield)



## USA & Canada

### LabLogic Systems, Inc.

East Pointe Park, 1040 East Brandon Blvd.  
Brandon, FL 33511-5509, USA

E-mail: [solutions@lablogic.com](mailto:solutions@lablogic.com)

Tel: +1-813-626-6848

Fax: +1-813-620-3708

Web: [www.lablogic.com](http://www.lablogic.com)



Certificate No: 10926  
ISO 9001

## Europe & Worldwide

### LabLogic Systems Limited

Paradigm House, 3 Melbourne Avenue  
Broomhill, Sheffield, S10 2QJ, UK

E-mail: [solutions@lablogic.com](mailto:solutions@lablogic.com)

Tel: +44 (0)114 266 7267

Fax: +44 (0)114 266 3944

Web: [www.lablogic.com](http://www.lablogic.com)



Certificate No: 1535  
ISO 9001



EXPERIENCE & EXPERTISE