



Specification

Hidex 300 SL

100-240V AC and 24 V DC.

Features Model (425-201) Automatic TDCR Liquid Scintillation Counter

General Description

Hidex 300SL automatic TDCR liquid scintillation counter for detection of low level beta and alpha radiation, suitable especially well for sophisticated research applications, direct DPM counting without External radioactive standard source and for extremely low levels of alpha activity when equipped with alpha/beta separation electronics.

The counter is equipped with unique 3 PMT technology to support TDCR counting and modern MS Windows software for user-friendly operation. The design with small footprint and low weight, allows integration into small laboratories and makes the system truly transportable.

The instrument has min. 70 mm internal Pb shielding design for low background.

There is no external radioactive standard needed. Transportation and movement of standard instrument does not require licenses for handling radioactivity and decommissioning is easy and low cost.

The instrument is intended for professional laboratory research use by trained personnel

Isotopes

Detection of beta and alpha radiation by Liquid Scintillation Counting (LSC)

- beta-emitters on LSC mode up to 2,000 keV.
- alpha-emitters on LSC mode up to 10,000 keV
- Preset protocols: H-3 (unquenched, high water load and low water load), C-14, S-35, P-33, P-32, I-125 and Free. Other nuclides and Cerenkov counting can be easily defined.
- Up to 4 energy regions (ROI) can be defined in same measurement

Analysis of data

- Multi-Channel Analysis
- 2 x logarithmic MCA with 1024 channels each
- Beta Spectrum Analysis using Report Render analysis program
- Quench curve analysis on MikroWin 2000 data reduction software or on MS Excel macro.
- Quench correction by TDCR method, External Eu-152 standard (option), quench curve method

Counting modes

- Beta (all coincidences), beta triple (triple coincidences), and alpha counting (a/b-separation, option)
- Direct measurement of beta efficiency by triple to double coincidence ratio
- Triple counting mode for luminescence free counting and for low level counting

Physical Dimensions

- 520 mm (Width) x 630 mm (Depth) x 680 mm (Height)
- Weight 125 kg

Operating conditions

- +10°C - +40°C; maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.
- Indoor use only
- Dust during operation: The system shall function in normal laboratory environment.
- Light during operation: The system shall function normally in normal ambient laboratory illumination

Storage and transport conditions

- Transportation conditions -40°C to +70°C, humidity 5 - 90 %, packed in transport packaging
- Storage conditions -25°C to +50°C, humidity 5 - 90 %, packed in transport packaging

Electrical Connections



Specifications

Hidex 300 SL automatic TDCR liquid scintillation counter

- Voltage: 100 ~ 240 VAC
- Frequency: 50 - 60 Hz
- Normal Power Consumption: <100W
- Max. Power Consumption: 350W

Samples

- Vials: 5 mL, 7 mL and 20 mL
- Vials require cap with flat surface to enable pneumatic gripping.
- Vial compatibility must be confirmed with Hidex.
- Vials are placed in a rack with 96 x 5 mL, 96 x 7ml vials and 40 x 20ml vials
- rack system is integrateable with liquid handlers and cell harvesters

Counting time

- 1 second - 168 h (7 d) /sample.
- Max number of repeats in one measurement protocol is 1024 (facilitates kinetic measurements)

Vial transport

- Robotic loading arm with pneumatic mechanism removes the need for a complex elevator mechanism.

Static electricity

- By deionizer prior the movement of the vial to detector.
- Deionization time can be defined by user at 1 second intervals

Quench correction by

The instrument has automatic colour, chemical and luminescence (quenching) detection and correction. Quench correction can be done using:

1. TDCR method

- Allows automatic direct DPM counting of single beta emitters without the use of External standard source or quench curve method
- Allows automatic direct DPM counting of samples labeled with dual beta emitters, and with constant quench, without the use of external radioactive standard source or quench curve method
- Deviation from true DPM is +/- 10 % (+/- 5 % typical) depending on the isotope and on the degree of quenching
- Accuracy of TDCR method can be increase down to +/- 2 % by combining with quench curve method

2. External Eu-152 standard (option)

- suggested only for double and triple labeled samples with variable quench
- requires fitting of the results on standard curve

3. Quench curve method

- set of quench standards required for each isotope measured
- can be used to improve the accuracy of TDCR method

4. Instant DPM

- by TDCR method

Instrument operation and software

Counter can be operated using:

- MS Windows based, CFR 21 part 11 compliant MikroWin 2000 user interface software,
- CommFiler CSV service software

Data reduction by:

- MS Windows based, CFR 21 part 11 compliant MikroWin 2000 software
- Excel based Report Render spectral analysis software
- By user definable export format (txt, xls, csv, etc.) to any other data reduction software or data base.

MikroWin 2000 parameter-files for measurement for most typical isotopes and applications are available at parameter library included in all instrument deliveries.

The parameter files include predefined settings for running the assay and calculation formulas for printing out desired data. Mikrowin 2000 software includes capabilities for biological assays, immunoassays, screening type assays, data validation and control history.



Specifications

Hidex 300 SL automatic TDCR liquid scintillation counter

Analysis of data includes features such as automatic calculation of activity, normalization of results, background subtraction, precision, error of measurement, minimum detectable activity, analysis of kinetic data, analysis of dual- and triple label data, half-life correction, etc. Also sample identification with sample number, position in the rack and with user defineable ID is facilitated.

Report Render spectral analysis software enables convenient and easy-to-use tools for optimizing FOM, analysis of spectra and background spectral subtraction. It also incorporates useful two-dimensional alpha-beta separation graph for validating alpha-beta separation performance.

MikroWin 2000 software is compatible with Windows 98, 2000, XP, Vista and Windows 7

Results are saved permanently in the data-file in computer memory after counting of each sample.

Communication connections

- Output port using the industry-standard RS-232 serial interface allows communication with PC.
- USB optional

Detectors

Three low background PMTs provide optimum measurement geometry and facilitate TDCR counting

- 3 x 1.5 inch Single Photon Counting Photo Multiplier Tube (PMT)
- Spectral range from 300 nm to 630 nm
- PMTs 120 degrees apart from each other

Detector shielding

- Optimal Pb shield design with a 4-pi 70 mm shielding of PMTs in all directions provides low background and minimizes instrument weight.
- Lead shutter provides optimum shielding from cosmic radiation
- Copper shielding around the measurement chamber eliminates x-rays from Pb shield.
- Measurement chamber with high reflective opaque paint maximizes light collection.

Performance Specifications

All the measurements are performed at temperature of $22^{\circ}\text{C} \pm 2^{\circ}\text{C}$ and at normal humidity conditions of Hidex laboratory, Turku, Finland (relative humidity not measured).

Counting efficiency:

- up to 73 % for H-3 and up to 98 % for C-14 with unquenched samples.
- > 30 % for H-3 quenched (8 mL water sample + 12 mL AquaLight cocktail)
- α 's (Po-210, U-234/238, Am-241, Rn-222, Ra-226) > 95 %

Typical background:

- < 12.5 CPM for H-3 in 8 mL water + 12 mL AquaLight (25 % eff.)
- 8.5 CPM for H-3 in 8 mL water + 12 mL AquaLight, with LL PMT option (25 % eff.)

FOM (E^2/B):

- H-3 in water 80 (> 90 with Low Level PMT option)
- H-3 unquenched, optimized window 180 (350 with Low Level PMT option)
- C-14 unquenched, optimized window 420 (800 with Low Level PMT option)

Count Rate

- up to 3.999.999 CPM, with MCA

Minimum detectable activity (MDA)

- Depends on factors such as the isotope, sample matrix, counting time and sigma value
- Lowest limit of detection for H-3 in water:
 - 17 Bq/l, 1 h counting time, standard unit/ 13 Bq/l with Low Level PMT option
 - 12.5 Bq/l, 2 h counting time, standard unit/ 9.2 Bq/l with Low Level PMT option

Performance Assessment

- Using standard samples
- OQ parameterfile with history log
- Standards with OQ history log can be included in all measurements



Specifications

Hidex 300 SL automatic TDCR liquid scintillation counter

Data recovery

Automatic storing of measured data in data-file created before the start of the measurement. The data of every sample is stored permanently in computer hard disk after completing the measurement of that particular sample.

Vial recovery

Automatic return of sample vials after power failure.

Power up diagnostics

- Automatic power up diagnostics of instrument settings.

Safety Standards

- CE marked
- External power UL/CSA approved

Warranty

- 12 months

Maintenance

- No recommended spares during the warranty period
- Extended warranty and maintenance & service program available on request.

21 CFR Part11 compliance

- Included with MikroWin 2000 software delivery

Equipment included in the delivery

- 425-201 Hidex 300 SL
- Owner's Handbook (on CD)
- Rack for 40 x 20 mL vials
- Rack for 96 x 7 mL vials
- Power supply 100-240V/24V
- RS-232 cable
- CommFilerCSV service software
- Excel based Report Render spectral analysis software
- MikroWin 2000 counter control and data reduction software with parameter file library for most common isotopes and applications (21 CFR part11 compliant)
- Parameter library for MikroWin software
- QC report

Optional Items

425-2001 Incubation specifications

- Cools down and stabilize the instrument 5°C below the ambient temperature

525-003 Alpha-beta separation

- High performance alpha beta separation.
- Easy to use two dimensional alpha beta separation plot for convenient separation optimization.
- No need to use pure standards to optimise separation conditions before counting
- Ultra low misclassification at 0.5 % on 100 000 CPM of activity
- Alpha background typically 0.3 CPM for Gross Alpha
- Alpha background typically < 0.1 CPM for specific energies

425-019 External Eu-152 standard

- External standard is needed with double label variable quench applications only.
- Low activity 74 kBq Eu-152. Under the IAEA free limit

425-018 Low Level PMTs

Counting efficiency:

- Up to 71 % for H-3, unquenched
- Up to 97 % for C-14, unquenched



Specifications

Hidex 300 SL automatic TDCR liquid scintillation counter

Typical background:

- 7.5 CPM for H-3 in 8 mL water + 12 mL AquaLight (25 eff.)
- 5 CPM for unquenched H-3 in optimized window
- 5 CPM for unquenched C-14 in optimized window

FOM (E^2/B):

- H-3 in water 90
- H-3 unquenched, optimized window 350
- C-14 unquenched, optimized window 800

PC, keyboard and printer

- Can be obtained locally
- Requirements of the PC:
 - MikroWin Hidex 2000 runs under Microsoft Windows 98, Windows NT 4.0, Windows 2000, Windows XP, Windows Vista and Windows 7 operating systems.
 - The computer requires a minimum of 32 MB available memory (64 MB or more is preferable).
 - Super VGA or higher resolution monitor is required (XGA is recommended)
 - Depending on installation options, a MikroWin Hidex 2000 installation requires 10 to 15 MB of available hard-disk space.
 - If there is not serial port in the PC, a USB to serial adapter must be included in the delivery to facilitate communication.
 - MS Excel is suggested to facilitate installation of Report Render spectral analysis software
 - Internet connection is suggested to facilitate updates from software manufacturer and remote desktop support from Hidex.

Manufacturer

Hidex Oy, Mustionkatu 2, FIN-20750 Turku, Finland

Data and specifications are subject to change, Hidex reserves the right to alter specifications.

No maximum is specified for background values.

Excel is trademark of Microsoft Corp.