COCKTAILS FOR LSC

Hidex provides high performance scintillation cocktails tested to give optimum performance with Hidex scintillation counters Hidex 300SL, Hidex 600SL and Hidex Triathler for aqueous, lipophilic, filter counting, alpha beta separation, general beta counting and low level applications. NPE free liquid scintillation cocktails are also available.

Hidex MaxiLight+

Hidex MaxiLight+ cocktail is a lipophilic, unquenched, DIN based cocktail with highest counting efficiency for organic and non-aqueous samples like bio-oils, dry samples, dried filters and for extraction of Rn-222 from water. MaxiLight+ has best alpha-beta separation efficiency. MaxiLight+ cannot be used for water samples or samples with any moisture.

Properties

- High performance cocktail for organic samples, biofuels and dried filters
- Excellent α/β -counting properties
- Optimum for Rn-222 counting with extractive method
- Does not accept any water or moisture
- NPE free
- Safe
- High flash point
- Odourless
- Low permeability through plastic vials

MaxiLight+

High performance cocktail for organic samples, bio-fuels and counting of dried filters. Excellent α/β -counting properties and optimum for Rn-222 counting with extractive method. Does not accept any water or moisture.

Product information: Specifically designed for high efficiency sample counting of organic samples. Safe, odourless, high counting efficiency, low permeation through plastic vials.

Contains: No NPE's (Nonylphenol Ethoxylates) EU WFD 2000/60/EC.

KEEP TIGHTLY CLOSED. PROTECT FROM SUNLIGHT. STORE AT OR SLIGHTLY BELOW ROOM TEMPERATURE.

ANGER

CATALOGUE NO.
CONTENTS (LITER):

461-003 461-015

1 5

FOR LABORATORY RESEARCH USE ONLY

H316 Causes skin irritation. H318 Causes serious eye damage. P264 Wash ... thoroughly after handling. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302-P352 IF ON SKIN: Wash with plenty of water. P305-P351-P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337-P313 If eye irritation persists: Get medical advice/attention.

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Storage

Store at or slightly below room temperature (15-25°C) in chemical storage cabinet. Keep container tightly sealed. Protect from light.

Transport information

UN3082, ENVIRONMENTALLY HAZARDOUS SUBSTANCE LIQUID, N.O.S (Di-isopropyl naphthalene isomers)
(IATA DGR: NOT RESTRICTED AS PER SPECIAL PROVISION A197)

Hidex AquaLight+

Hidex AquaLight+ is a general-purpose high safety scintillation cocktail for aqueous and non-aqueous samples, for both general beta counting and alpha-beta separation applications. AquaLight+ is a DIN based cocktail capable of handling a broad range of solutes, combining high counting efficiency and low background with a very high level of sample acceptance, particularly for high ionic strength solutes.

Properties

- Water-dissolving DIN cocktail for general beta counting and α/β separation applications
- High sample capacity, low background and high counting efficiency
- Solubility up to 12 ml water in 10 ml cocktail
- Compatible with urine samples and samples with high mineral acid concentrations.
- Good α/β -counting properties
- No diffusion through plastic vials ensures stable repeatable counting
- High flash point of ~140°C increases safety
- REACH regulations classify NPEs used in this cocktail as SVHC's (Substances of Very High Concern) and therefore this cocktail is available for research use only.
- Disposal by incineration only

General purpose cocktail for aqueous samples with high sample capacity, low background and high counting efficiency. Compatible with urine samples and samples with high mineral acid concentrations. Good α/β -counting properties.

Product information: Low background, high counting efficiency. Colour & quench resistant. Low permeation through plastic vials.

KEEP TIGHTLY CLOSED. PROTECT FROM SUNLIGHT. STORE AT OR SLIGHTLY BELOW ROOM TEMPERATURE.

DANGER

CATALOGUE NO.	CONTENTS (LITER):
461-004	1
461-006	5

FOR LABORATORY RESEARCH USE ONLY

EN: H504 May be fatal if swallowed and enters airways. H302 Causes skin irritation. H331 Causes eye irritation. H403 May cause long-term harmful effects to aquatic life. P201 Do not breathe vapour. P273 Do not release to the environment. P280 Wear eye protection, gloves, protective clothing and eye face protection. P301+P312 IF SWALLOWED, immediately call a POISON CENTRE or doctor/physician. P302+P352 IF ON SKIN Wash with plenty of water. P303+P361+P353 IF IN EYES, rinse copiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing. P304+P340 If you feel unwell, get medical attention. P305+P351+P338 IF eye irritation persists; Get medical advice/attention. P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

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Storage

Store at or slightly below room temperature (15-25°C) in chemical storage cabinet. Keep container tightly sealed. Protect from light.

Transport information

UN3082, 3082 ENVIRONMENTALLY HAZARDOUS SUBSTANCE LIQUID, N.O.S. (Contains Di-isopropyl-naphthalene isomers, 4-Nonyl phenol, branched ethoxylated)
(IATA DGR: NOT RESTRICTED AS PER SPECIAL PROVISION A197)

Hidex AquaLight+ Ultra Low Level

Hidex AquaLight+ Ultra Low Level cocktail is water-dissolving DIN based cocktail with highest performance for low level tritium counting, high sample capacity, ultra low background and high efficiency at low temperatures.

Properties

- Water-dissolving high performance cocktail for low level tritium counting
- High sample capacity, ultra low background and high counting efficiency at low temperatures
- Solubility up to 12 ml water in 10 ml cocktail
- Stable at temperatures down to 12°C for sub-ambient counting
- Compatible with urine samples and samples with high mineral acid concentrations.
- Good α/β -counting properties
- No diffusion through plastic vials ensures stable repeatable counting
- High flash point of ~140°C increases safety
- REACH regulations classify NPEs used in this cocktail as SVHC's (Substances of Very High Concern) and therefore this cocktail is available for research use only.
- Disposal by incineration only

AquaLight+ Ultra Low Level

FOR LABORATORY RESEARCH USE ONLY

High performance cocktail for low level tritium counting with high sample capacity, ultra low background and high counting efficiency at low temperatures.

Product information: Low background, high counting efficiency. Colour & quench resistant. Low permeation through plastic vials

KEEP TIGHTLY CLOSED. PROTECT FROM SUNLIGHT. STORE AT OR SLIGHTLY BELOW ROOM TEMPERATURE.

DANGER

CATALOGUE NO.
 CONTENTS (LITER):

461-007
5

H304 May be fatal if swallowed and enters airways. H315 Causes skin irritation. H330 Causes eye irritation. H410 May cause long-lasting harmful effects to aquatic life. P201 Do not breathe vapour. P273 Avoid release to the environment. P280 Wear suitable protective clothing, gloves and eye/face protection. P301-P310 IF SWALLOWED: Immediately call a POISON CENTRE or doctor/physician. P302-P352 IF ON SKIN: Wash with plenty of water. P303-P361-P353 IF IN EYES: rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing. P304 Do not induce vomiting. P305-P351-P338 In case of eye irritation: lift medical advice/attention. P307-P311 In case of eye irritation persists: lift medical advice/attention. P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

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Cat. No. 461-007 AquaLight+ Ultra Low Level 5 L plastic bottle

Storage

Store at or slightly below room temperature (15-25°C) in chemical storage cabinet. Keep container tightly sealed. Protect from light.

Transport information

UN3082, 3082 ENVIRONMENTALLY HAZARDOUS SUBSTANCE LIQUID, N.O.S. (Contains Di-isopropylNaphthalene isomers, 4-Nonyl phenol, branched ethoxylated)
 (IATA DGR: NOT RESTRICTED AS PER SPECIAL PROVISION A197)

Hidex AquaLight AB

Hidex AquaLight AB is a water-dissolving DIN based cocktail for alpha-beta separation applications, such as direct detection of Rn-222 in water or gross alpha-beta measurements with good alpha-beta separation efficiency. Solubility up to 6.5 mL water in 13.5 mL cocktail.

Properties

- Water-dissolving DIN cocktail for α/β separation applications
- High alpha-beta separation efficiency
- Solubility up to 6.5 ml water in 13.5 ml cocktail
- Good α/β -counting properties
- No diffusion through plastic vials ensures stable repeatable counting
- No REACH regulations classify NPEs used in this cocktail as SVHC's (Substances of Very High Concern) therefore this cocktail is available for research use only.
- Disposal by incineration only

AquaLight AB is high performance DIN based LSC cocktail for aqueous and non-aqueous samples. Accepts up to 6.5 ml water in 13.5 ml cocktail at 20 °C Optimized for α/β separation.

Product information: Odourless, colour and quench resistant. Low photo and chemiluminescence.

Contains: Diisopropyl naphthalene isomers.

KEEP TIGHTLY CLOSED. PROTECT FROM SUNLIGHT. STORE AT OR SLIGHTLY BELOW ROOM TEMPERATURE.

ANGER

CATALOGUE NO.
CONTENTS (LITER):

461-021 461-020

1 **5**

FOR LABORATORY RESEARCH USE ONLY

H304 May be fatal if swallowed and enters airways. H315 Causes skin irritation. H319 Causes serious eye irritation. H400 Very toxic to aquatic life with long lasting effects. H410 Very toxic to aquatic life with long lasting effects. P273 Avoid release into the environment. P280 Wear protective gloves/protective clothing/eye protection/face protection. P301 + P310 If SWALLOWED: Immediately call a POISON CENTRE / DOCTOR. P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. P501 Dispose of contents / container in accordance with local/national regulations.

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Storage

Store at or slightly below room temperature (15-25°C) in chemical storage cabinet. Keep container tightly sealed. Protect from light.

Transport information

UN3082, 3082 ENVIRONMENTALLY HAZARDOUS SUBSTANCE LIQUID, N.O.S. (Contains Di-isopropyl naphthalene isomers, 4-Nonyl phenol, branched ethoxylated)
(IATA DGR: NOT RESTRICTED AS PER SPECIAL PROVISION A197)

Hidex AquaLight Beta

AquaLight Beta is DIN based multi-purpose LSC cocktail for aqueous and nonaqueous samples. Specifically designed for high sample capacity with high efficiency.

Properties

- Optimized for general β -counting
- Accepts up to 10 ml of dilute aqueous samples in 10 ml cocktail at 20 °C
- Extremely resistant to chemiluminescence and is therefore suitable for use with the solubilised samples.
- REACH regulations classify NPEs used in this cocktail as SVHC's (Substances of Very High Concern) and therefore this cocktail is only available for research use only.

AquaLight Beta is DIN based LSC cocktail for aqueous and non-aqueous samples. Specifically designed for high sample capacity with high efficiency. Accepts up to 10 ml water in 10 ml cocktail at 20 °C Optimized for β -counting.

Product information: Odourless, colour and quench resistant.

Contains: Diisopropyl naphthalene isomers.

KEEP TIGHTLY CLOSED. PROTECT FROM SUNLIGHT. STORE AT OR SLIGHTLY BELOW ROOM TEMPERATURE.

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Cat No. 461-035 AquaLight Beta 5 L plastic bottle

Storage

Store at or slightly below room temperature (15-25°C) in chemical storage cabinet. Keep container tightly sealed. Protect from light.

Transport information

UN3082, 3082 ENVIRONMENTALLY HAZARDOUS SUBSTANCE LIQUID, N.O.S. (Contains Di-isopropyl naphthalene isomers, 4-Nonyl phenol, branched ethoxylated)
(IATA DGR: NOT RESTRICTED AS PER SPECIAL PROVISION A197)

Oxidizer 600 OX Radiocarbon

Oxidizer 600 OX Radiocarbon is a complete liquid scintillation cocktail for trapping and counting radioactive carbon dioxide. This cocktail is designed for the Hidex 600 OX Oxidizer. The cocktail contains a carbon dioxide absorber that effectively traps CO₂ into the cocktail. The cocktail contains also scintillator and, therefore, the combusted sample in the cocktail is instantly ready for the LSC measurement.

Properties

- 1 ml of cocktail absorbs 2.8 mmol of carbon dioxide
- High counting efficiency and low background
- Effective carbon dioxide absorption by 3-methoxypropylamine
- Contains scintillator for instant C-14 measurement
- NPE free
- Disposal by incineration only

Oxidizer 600 OX Radiocarbon

A complete liquid scintillation cocktail for trapping and counting radioactive carbon dioxide. To be used with Hidex 600 OX Oxidizer.

Product information: 1 ml of cocktail absorbs 2.8 mmol of carbon dioxide.

KEEP TIGHTLY CLOSED. PROTECT FROM SUNLIGHT. STORE AT OR SLIGHTLY BELOW ROOM TEMPERATURE.

CATALOGUE NO.

461-071

FOR LABORATORY RESEARCH USE ONLY

H225 Flammable liquid and vapour. H232 Harmful if inhaled. H315 Causes skin irritation. H318 Causes serious eye irritation. H335 H336 May cause respiratory irritation. May cause drowsiness or dizziness. H411 Toxic to aquatic life with long lasting effects. P201 Keep away from heat/sparks/open flames/hot surfaces – No smoking. P201 P202 Wear protective gloves/protective clothing/eye protection/face protection. P203–P211–P231 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. P203–P211–P231 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. P204–P240 IF INHALED: Remove person to fresh air and keep comfortable for breathing.

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Ordering Info

Cat. No. 461-071 Oxidizer 600 OX Radiocarbon 1 L

Delivered in glass bottles, bottle size depends on order volume and transport method.

Storage

Store at or slightly below room temperature (15–25°C) in chemical storage cabinet. Keep container tightly sealed. Protect from light.

Transport information

UN2924 FLAMMABLE LIQUID, CORROSIVE, N.O.S. (3-methoxypropylamine, 1,2,4-trimethylbenzene, isobutyl alcohol) ENVIRONMENTALLY HAZARDOUS)

Oxidizer 600 OX Tritium

Oxidizer 600 OX Tritium is a complete liquid scintillation cocktail for trapping and counting tritium-water from combustion systems. This cocktail is an NPE free pseudocumene based LSC cocktail especially for the Hidex 600 OX Oxidizer. The cocktail is designed to accept water from sample oxidisers without foaming or gelling and produces a stable microemulsion suitable for LSC counting.

Properties

- High volume water capacity, accepts up to 6 ml of water into 10 ml of cocktail
- High tritium counting efficiency
- No foaming
- NPE free
- Disposal by incineration only

A liquid scintillation cocktail for counting tritium produced from sample oxidation. To be used with Hidex 600 OX Oxidizer.

Product information: Accepts up to 6 ml of water into 10 ml of cocktail. NPE free.

Storage information: Keep tightly closed. Protect from sunlight. Store at or slightly below room temperature. Keep away from heat.

KEEP AWAY FROM HEAT
KEEP TIGHTLY CLOSED. PROTECT FROM SUNLIGHT. STORE AT OR SLIGHTLY BELOW ROOM TEMPERATURE.

CATALOGUE NO.

461-073

FOR LABORATORY RESEARCH USE ONLY

H225 Flammable liquid and vapour. H316 Causes skin irritation. H318 Causes serious eye damage. H332-H361-H411 Harmful if inhaled. May cause respiratory irritation. Toxic to aquatic life with long lasting effects. P201 P202 Keep away from heat/sparks/open flames/hot surfaces – No smoking. P273 Avoid release to the environment. P280 Wear protective gloves/protective clothing/eye protection/face protection. P303-P361-P353 IF ON SKIN (or hair): Remove/ take off immediately all contaminated clothing. Rinse skin with water/shower. P304-P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing. P305-P351-P338 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing.

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Ordering Info

Cat. No. 461-073 Oxidizer 600 OX Tritium 1 L

Delivered in glass bottles, bottle size depends on order volume and transport method.

Storage

Store at or slightly below room temperature (15-25°C) in chemical storage cabinet. Keep container tightly sealed. Protect from light.

Transport information

UN1993 FLAMMABLE LIQUID, N.O.S. (1,2,4-Trimethylbenzene), ENVIRONMENTALLY HAZARDOUS

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Manufactured by Acton Ltd • 03100012010 • Contact: 01473 841100