

Laura RadiopharmaTM

PET/SPECT Radiochromatography Instrument Control,
Data Collection and Analysis Software

www.lablogic.com



LabLogic
EXPERIENCE & EXPERTISE

The ideal solution for PET/SPECT facilities

Laura Radiopharma™ is the industry standard single point control and data analysis software. It supports over 350 r-HPLC detectors, r-TLC scanners, MCA, GC and HPLC modules from leading manufacturers.

Whether you work within a regulated or academic environment, you will only need to learn one software system for the complete suite of PET/SPECT QC instruments.

Easy to use

Laura Radiopharma™ uses recognisable formats and wizards for multi-function tasks to make it simple to use by multiple users with different levels of experience with the system.

A single software solution for the range of disciplines needed in the QC department makes the learning process easy and the day to day use effective.

Compliant

Laura Radiopharma™ can be installed and configured to meet the most rigorous regulatory compliance requirements. Featuring configurable audit trail, multi level security, e-signatures and data storage into a secure database environment Laura Radiopharma™ is the PET/SPECT industry's answer for regulatory compliance.

Connectivity

Laura Radiopharma™ provides connectivity to the widest range of detectors and chromatography systems, including the Flow-RAM and Scan-RAM for radioactivity detection along with chromatography systems for HPLC, GC as well as MCA.

Prudent Investment

From standalone installation to a full client server implementation, Laura Radiopharma™ offers a scalable solution to your data storage requirements. Application and data can be separated across local clients and central servers to match IT requirements.

Consistent and Secure


The software allows the user to lock methods, reports, etc., to ensure consistency of data.

Digital Signal

Signal from radiochromatography instrumentation is reported digitally in counts, CPS and CPM.


Half-Life Correction

Half-life correction function enables users to correct for half life while a run is in progress or post run using a reference time and date.



Reporting

Relevant values reported such as % total for radiochemical purity checks. The powerful Report Designer allows you to configure multiple reports, calculate and present data in real time.



Regulatory Compliance


Laura Radiopharma™ is designed to meet the regulatory compliance requirements for GxP and FDA 21 CFR Part 11 and FDA 21 CFR Part 212 requirements.

Complete Workflow Solution

Laura Radiopharma™ will guide the analyst through the complete QC procedure. Methods can be pre-defined and chosen from a drop down list. Automatic settings provide for automatic peak integration and report printing.

Intuitive

Laura Radiopharma™ is easy to use allowing quick familiarisation with the system and rapid adoption by its users.



Single Software Solution

Laura Radiopharma™ is a single software solution for the PET/SPECT QC environment. Rather than having several software systems to use, the QC analyst can use Laura Radiopharma™ for r-HPLC, HPLC, r-TLC, GC and MCA analysis.

Service and Support

Users of our systems can benefit from our comprehensive, fully inclusive service and support.

We can give reassurance that if things go wrong or you need expert advice, help is only an e-mail or phone call away.



Validation Services

Our Validation Service enables you to implement and get maximum value from your investments as soon as possible.

We work as a partner with your Quality Manager, System Manager and users to provide a tailored Validation Plan, suited to your needs. Our Validation Specialists who have many years' of experience in GLP system validation, detailed knowledge of our systems, together with other industry standard systems to help you meet company and regulatory requirements.

Training

LabLogic can provide a variety of training courses and workshops to help you get the most out of your instrument and software.

All training is performed by our expert Product and Support Specialists who have many years experience in the development and use of the instruments and software.

Certificates can be provided to complement your internal GLP training records.

USA & Canada

LabLogic Systems, Inc.
3901 Centerview Drive, Suite B
Chantilly, VA 20151, USA

E-mail: solutions@lablogic.com
Tel: +1-703-429-4209

www.lablogic.com

Europe & Worldwide

LabLogic Systems Limited
Innovation House, 6 Europa View
Sheffield, S9 1XH, UK

E-mail: solutions@lablogic.com
Tel: +44 (0)114 266 7267

www.lablogic.com



Certificate No: 10926
ISO 9001



Certificate No: 1536
ISO 9001

INVESTORS IN PEOPLE
We invest in people Gold