



Laura for PETTM

PET/SPECT Radiochromatography
Instrument Control, Data Collection
and Analysis Software

www.lablogic.com



LabLogic
EXPERIENCE & EXPERTISE

The ideal solution for PET/SPECT facilities

Laura for PET is the industry standard single point control and data analysis software. It supports over 350 r-HPLC detectors, r-TLC scanners, MCA, GC and HPLC modules from leading manufacturers.

Whether you work within a regulated or academic environment, you will only need to learn one software system for the complete suite of PET/SPECT QC instruments.

Easy to use

Laura for PET uses recognisable formats and wizards for multi-function tasks to make it simple to use by multiple users with different levels of experience with the system.


A single software solution for the range of disciplines needed in the QC department makes the learning process easy and the day to day use effective.

Compliant

Laura for PET can be installed and configured to meet the most rigorous regulatory compliance requirements. Featuring configurable audit trail, multi level security, e-signatures and data storage into a secure database environment Laura for PET is the PET/SPECT industry's answer for regulatory compliance.

Connectivity

Laura for PET provides connectivity to the widest range of detectors and chromatography systems, including the Flow-RAM and Scan-RAM for radioactivity detection along with chromatography systems for HPLC, GC as well as MCA.



Service and Support

Users of our systems can benefit from our comprehensive, fully inclusive service and support.

We can give reassurance that if things go wrong or you need expert advice, help is only an e-mail or phone call away.

Validation Services

Our Validation Service enables you to implement and get maximum value from your investments as soon as possible.

We work as a partner with your Quality Manager, System Manager and users to provide a tailored Validation Plan, suited to your needs. Our Validation Specialists have years of experience in GLP system validation, detailed knowledge of our systems, together with other industry standard systems to help you meet company and regulatory requirements.

Training

LabLogic can provide a variety of training courses and workshops to help you get the most out of your instrument and software.

All training is performed by our expert Product and Support Specialists who have many years experience in the development and use of the instruments and software.

Certificates can be provided to complement your internal GLP training records.

USA & Canada

LabLogic Systems, Inc.
1911 N US HWY 301, Suite 140
Tampa, FL 33619, USA

E-mail: solutions@lablogic.com
Tel: +1-813-626-6848
Fax: +1-813-620-3708

www.lablogic.com

Europe & Worldwide

LabLogic Systems Limited
Paradigm House, 3 Melbourne Avenue
Broomhill, Sheffield, S10 2QJ, UK

E-mail: solutions@lablogic.com
Tel: +44 (0)114 266 7267
Fax: +44 (0)114 266 3944

www.lablogic.com

