

Product Datasheet

Radhound X/I

Hand-held Digital Radiation Meter

The Radhound X/I is an advanced hand-held survey meter with an internal dose rate detector for x-ray and gamma measurement.



Application Areas

- Health physics: for contamination monitoring on surfaces, clothing and objects etc.
- Nuclear medicine departments: suitable for ^{125}I , $^{99\text{m}}\text{Tc}$, etc.
- Radiological survey work and lab use.
- Emergency planning, response and clean up.
- Research applications.

Features

- Clear digital LCD display with backlight.
- GM detector for contamination and dose rate measurements.
- Fully adjustable alarm levels.
- Scaler timer function.
- Peak mode.
- Over range.
- Dose rate range: 0.01 $\mu\text{Sv/hr}$ – 1 mSv/hr

Specifications

Mechanical

Dimensions 165 x 115 x 59 mm (excluding handle)

Weight Typically 1.4 kg

Units uSv/hr, mSv/hr, Sv/hr

Display Large backlit LCD 70 mm x 40 mm
Automatic backlight movement sensor – when instrument is picked up, the backlight comes on
On screen battery status indication

Main Modes of Operation Dose rate, integration mode, histogram, peak value
Menu gives control of system parameters and multiple probe library

Radiological Performance

Measured Quantity Ambient Dose equivalent H*(10)

Dose Rate Range 0.01 uSv/hr – 1 mSv/hr

Energy Response 20 keV – 1.5 MeV

Power

Battery Uses 2 x 1.5 V Alkaline C-Cell (or optional Lithium Ion rechargeable)

Battery Life Typically >8 hours continuous use

Environmental

Operating Temperature -10°C to +50°C

Storage Temperature -25°C to +60°C

Compliance

Standard Compliance In conformity with EMC directive (2004/108/EC), EN61326-1
EN61000-3-2, EN6100-3-3
Low Voltage Directive (2006/95/EC)
Designed to meet IEC 60325-2006 and IEC 60846-2004 Nuclear Instrumentation Standards

Specifications are subject to change without notice.
For the most up-to-date specifications, please visit www.southernscientific.co.uk

Europe & Worldwide

LabLogic Systems Limited
Innovation House, 6 Europa View
Sheffield, S9 1XH, UK

E-mail: solutions@lablogic.com

Tel: +44 (0)114 266 7267

www.lablogic.com

USA & Canada

LabLogic Systems, Inc.
3901 Centerview Drive, Suite B
Chantilly, VA 20151, USA

E-mail: solutions@lablogic.com

Tel: +1-703-429-4209

www.lablogic.com



INVESTORS IN PEOPLE
We invest in people Gold

