

^{18}O Water

Patented ultra-pure oxygen distilled ^{18}O water



Superior Quality

Stable Supply

Consistency

Manufactured by Taiyo Nippon Sanso in Japan, our ^{18}O Water is produced under cGMP conditions using a patented ultra-pure oxygen distillation process. This unique process results in a superior quality target material with the highest enrichment available at >98atom% and a chemical purity of >99.99%. Typical benefits include greater yields from synthesis modules and higher specific activity.

Specifications

Formula	H_2^{18}O
Enrichment	>98atom%
Chemical Purity	>99.99%
Molecular Weight	19.98
Container	Glass vial with septum top and aluminum crimp seal
Package Size	10g, 20g, and 50g vials

LabLogic Systems, Inc.



Certificate No: 10926
ISO 9001

East Pointe Park, 1040 East Brandon Blvd.
Brandon, FL 33511-5509, USA
E-mail: solutions@lablogic.com
Tel: +1 (813) 626-6848
Web: www.lablogic.com



EXPERIENCE & EXPERTISE